

Detection and Quantification of Tert-Butyl Mercaptan (TBM) in a Natural Gas Matrix Using a Transportable Micro Gas Chromatography System with a Thermal Conductivity Detector



Presented by Debbie Hutt

PITTCON Conference and Expo 2012

Tuesday, March 13, 2012

3:35PM

Advances in Fuels and Petrochemical Analysis II

Presentation Outline

- Odorant Background
- Why Monitor Odorants?
- Odorant Regulations
- TBM and THT as Odorants
- The Need for Precise Analysis
- Instrumentation
 - 3000 Micro GC Highlights
 - Modular Design
 - Thermal Conductivity Detector (TCD)
- Data and Repeatability




Odorant Background

- Natural gas is colorless and odorless
- Odorants are added to natural gas pipelines to serve as a leak indicator and warning agent
- Odorants are sold as blends of two or more mercaptans, alkyl sulfides, or cyclic sulfides
- Common odorants:
 - TBM – *Tert-Butyl Mercaptan*
 - IPM – Isopropyl Mercaptan
 - NPM – Normal Propyl Mercaptan
 - SBM – Secondary Butyl Mercaptan
 - THT – *Tetrahydrothiophene/Thiophane*
 - DMS – Dimethyl Sulfide
 - MES – Methyl Ethyl Sulfide



Why Monitor Odorants?


- Many factors contribute to “**Odor Fading**” which will reduce the amount or intensity of the odorant over time
- Factors include:
 - Soil adsorption
 - Pipeline adsorption/absorption from plastic or steel pipes
 - Contaminants in “wet” natural gas (high amount of C4+ HCs) can mask odorants
 - Oxidation (air or rust are present) of mercaptans to disulfides



- Presence of odorants such as ethyl and methyl mercaptan can convert added mercaptans (TBM, IPM, NPM) to disulfides
- Thus, concentrations of odorant must be monitored frequently and consistently


Odorant Regulations

- The Department of Transportation (DOT) requires odorants to be added in accordance with **Federal Regulation 192.625**
- There must be enough odorant so that a person with a “normal” sense of smell can detect it
- **192.625** does not specify which blend, how often to sample, or where to sample
 - It is left up to the operator to determine the blend, frequency, and location
 - Each natural gas operation is different
 - Factors that come into play:
 - Previous leak history
 - Age of pipelines
 - Presence of new lines in the system
 - Type of piping
 - “Wet” vs “Dry” natural gas
 - Naturally odorized vs. synthetically odorized




TBM as an Odorant

- TBM is one of the most common odorants in the US
 - Highest resistance to oxidation (out of the mercaptans)
 - Low odor threshold (<0.1 ppb)
 - Excellent soil penetration
 - **High freezing point (34°F, 1.1°C), so it must be blended with other odorants**
 - A common blend is 80% TBM, 15% IPM, and 5% NPM
 - Industry expectation is ~1 ppm in the US



THT as an Odorant

- THT is a common odorant in Europe
 - Does not oxidize
 - **Low odor impact, not as pungent as TBM**
 - **Poor soil penetration**
 - Common blends are 100% THT, or 50% THT with 50% TBM
 - Industry expectation is **~4 ppm** in Europe



The Need for Precise Analysis

▪ Requirements:

- Separation and quantification of TBM and/or THT at **low ppm levels** in a natural gas matrix
- Simple and easy to use
- Transportable or portable instrumentation
- Precision



Gas Chromatography as an Analysis Option

- Gas chromatography (GC) technology provides separation and analysis capabilities for TBM and/or THT in a natural gas matrix
- GC software provides users with component concentrations
- Thermal conductivity detectors (TCDs) are universal detectors that are simple to operate and meet the required sensitivity needs

Instrumentation

- **INFICON 3000 Micro GC (MGC)**
 - MEMS based TCDs and injectors
 - 0.5 to 4 mL/min flow rates

- **1-Channel Configuration (for TBM)**
 - **Channel A: 15 m 14% Phenyl, 86% Dimethylpolysiloxane (CPSil13CB)**
 - 0.15 mm I.D.
 - Medium polarity
 - Large Variable Volume Injector (LVI) to provide the best sensitivity
 - Helium carrier gas

- **Additionally, a second channel can be utilized to analyze THT:**
 - 4 m 14% Cyanopropyl-phenyl, 86% Dimethylpolysiloxane (CPSil19CB, OV-1701)
 - LVI Injector



3000 Micro GC Highlights

- **User-friendly design**
 - Simple
 - Out-of-the-box ready
 - Rugged for field use


- **Size**
 - Small footprint
 - Lightweight
 - 2-Channel weighs 18 lbs
 - Low carrier gas consumption
 - < 5 cc/min per channel

3000 Micro GC Highlight, cont.

- **Fast**
 - Typical analysis times range from 1 to 3 minutes
 - TBM and THT can be analyzed in less than 90 seconds




- **Precise**
 - Precision equivalent to bench-top GCs

- **Transportable**
 - Easy to transport in the lab or in the field
 - On site sampling reduces the time and cost of sending samples to a lab



Modular Design Delivers Flexibility

- Each module is a self-contained GC comprising of a:
 - **Micro**-machined injector
 - High-resolution capillary column
 - **Micro**-machined TCD



- Modules can easily be replaced within minutes to adapt quickly to different application needs
- Minimal downtime

Chromatographic Module


Injector

Capillary Column

TCD



**Electronic "D board"
attached to back
of module**

**Digital
Pneumatics
(EPC)**




Thermal Conductivity Detector (TCD)

- Low internal volume
- Out-of-the-box ppm sensitivity
- Combined with LVI will give best sensitivity
- Dip switch for plug and play swapping option
- Linear from low ppm to 100%




TBM in Natural Gas

- Calibration gas:
 - 2ppm TBM in nitrogen
- The 3000 Micro GC was connected to a water heater supplied with **natural gas**
- Expected TBM levels are ~1 to 3 ppm



3000 Micro GC Method Parameters




Low temperature

High inject time

High detector sensitivity


TBM Repeatability in Natural Gas– 10 Runs

- TBM concentration was ~2.5 ppm




Channel	Number of Analyte	Compound	Retention Time	Area %RSD
A	1	TBM	62.44	2.581

TBM Repeatability in Natural Gas – 10 Run Overlay




TBM Limit of Detection

- 2 ppm TBM standard was diluted with a gas-tight syringe into a natural gas standard containing C6+ components
- ~0.75 ppm of TBM was analyzed



THT in Natural Gas

- Calibration gas:
 - 5 ppm THT in methane
 - 14% Cyanopropyl-phenyl, 86% Dimethylpolysiloxane column



Conclusion

- Using 1 or 2 channels, the 3000 Micro GC is an excellent instrument to analyze low ppm levels of TBM and THT odorants
- Within 90 seconds, the complete separation TBM and THT was achieved
- 2.5 ppm TBM has a repeatability of less than 3 %RSD over the course of 10 sequential runs
- The 3000 Micro GC was able to detect down to **0.75 ppm** of TBM in a natural gas matrix



Visit booth #3133

or

**Visit our website www.INFICON.com for more
information**

Questions?

